


Anti-Ass1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P00966

Item Number: CM0010943



Introduction

High-quality anti-ASS1 rabbit monoclonal antibody for WB, IF-P, and ELISA applications. Cross-reacts with human, mouse, and rat ASS1; assists in urea cycle and arginine metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ASS1 Rabbit mAb
Remarks/Alias	ASS; CTLN1; ASS1
Species	Human
Gene ID (Human)	445
Gene ID	445
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ASS1 (P00966).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00966
Protein Weight	47 kDa
Shipping	Ice bag