

Anti-Asun/Ints13 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nvm9

Item Number: CM0007394



Introduction

Anti-ASUN/INTS13 (Integrator complex subunit 13) Rabbit Monoclonal Antibody, validated for WB and ELISA. Recognizes human ASUN/INTS13, a key component of the integrator complex involved in transcription termination and cell cycle regulation. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	ASUN/INTS13 Rabbit mAb
Remarks/Alias	ASUN; GCT1; NET48; Mat89Bb; SPATA30; C12orf11; ASUN/INTS13
Species	Human
GeneID (Human)	55726
GeneID	55726
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human ASUN/INTS13. (NP_060634.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NVM9
Protein Weight	80kDa
Shipping	Ice bag