

Anti-Asymmetric Dimethyl-Histone H3-R17 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, If/Icc, Elisa

Item Number: CM0000282



Introduction

Rabbit polyclonal antibody targeting asymmetric dimethyl-histone H3 at arginine 17 (H3R17me2a). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat; predicted wide cross-reactivity. Core nucleosome component essential for epigenetic regulation.

[Learn More](#)

Parameter	Description
Product Name	Asymmetric DiMethyl-Histone H3-R17 Rabbit pAb
Remarks/Alias	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3Fj; H3C10; H3C11; HIST1H3J; Asymmetric DiMethyl-Histone H3-R17
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag