

Anti-At1R/Agtr1 Rabbit Polyclonal Antibody For Wb, Elisa - P30556

Item Number: CM0017152



Introduction

Rabbit polyclonal antibody targeting human AT1R/AGTR1, suitable for western blot and ELISA. Cross-reacts with mouse and rat. Recognizes the angiotensin II type-1 receptor, a key regulator of blood pressure. ~41 kDa.

[Learn More](#)

Parameter	Value
Product Name	AT1R/AGTR1 Rabbit pAb
Remarks/Alias	AT1; AG2S; AT1B; AT1R; ATR1; AT1AR; AT1BR; AT2R1; HAT1R; AGTR1B; AT1R/AGTR1
Species	Human
Gene ID (Human)	185
Gene ID	185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 294-388 of human AT1R/AGTR1 (NP_114438.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30556
Protein Weight	41kDa
Shipping	Ice bag