

Anti-Atad3A Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nvi7

Item Number: CM0029582

Introduction

Rabbit polyclonal antibody targeting human ATAD3A, suitable for WB, IHC-P, and ELISA applications. Validated for mitochondrial biology studies, including network organization and metabolism. Cross-reacts with mouse and rat. Ideal for cancer and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	ATAD3A Rabbit pAb
Remarks/Alias	ATAD3A; HAYOS; ATPase family; ATAD3A/B
Species	Human
Gene ID (Human)	55210
Gene ID	55210
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ATAD3A (NP_001164006.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVI7
Protein Weight	71kDa
Shipping	Ice bag