

Anti-Atad3A/Atad3B Rabbit Polyclonal Antibody For Wb, If/Icc, And Elisa - Q9Nvi7 / Q5T9A4

Item Number: CM0021642

Introduction

Rabbit polyclonal antibody targeting ATAD3A and ATAD3B, validated for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for mitochondrial research and cell metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	ATAD3A/ATAD3B Rabbit pAb
Remarks/Alias	TOB3; AAA-TOB3; ATAD3A/ATAD3B
Species	Human
GeneID (Human)	55210/83858
GeneID	55210/83858
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVI7/Q5T9A4
Protein Weight	71/73kDa
Shipping	Ice bag