

Anti-Atad3A/B Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nvi7 / Q5T9A4

Item Number: CM0006637

Introduction

Rabbit monoclonal antibody targeting ATAD3A/B, suitable for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Essential for mitochondrial network organization and metabolism; overexpressed in lung adenocarcinomas.

[Learn More](#)

Parameter	Value
Product Name	ATAD3A/B Rabbit mAb
Remarks/Alias	ATAD3A; HAYOS; ATPase family; ATAD3A/B
Species	Human
GeneID (Human)	55210/83858
GeneID	55210/83858
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ATAD3A/B. (NP_060658.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NVI7/Q5T9A4
Protein Weight	57kDa/66kDa/71kDa
Shipping	Ice bag