

Anti-Ataxin-3 (Atxn3) Rabbit Monoclonal Antibody For Wb And Elisa - P54252

Item Number: CM0006632



Introduction

Anti-Ataxin-3 (ATXN3) rabbit monoclonal antibody, validated for Western blotting (WB) and ELISA, with cross-reactivity to mouse and rat. Ideal for research on neurodegeneration, protein homeostasis, and deubiquitination.

[Learn More](#)

Parameter	Value
Product Name	Ataxin-3 (ATXN3) Rabbit mAb
Remarks/Alias	AT3; JOS; MJD; ATX3; MJD1; SCA3; Ataxin-3 (ATXN3)
Species	Human
GeneID (Human)	4287
GeneID	4287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-361 of human Ataxin-3 (ATXN3) (NP_004984.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P54252
Protein Weight	41kDa
Shipping	Ice bag