

Anti-Ataxin-3 (Atxn3) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P54252

Item Number: CM0014644



Introduction

Rabbit polyclonal antibody against human Ataxin-3, validated for WB, IF/ICC, ELISA. Reacts with human and mouse ATXN3, targeting amino acids 1-361. Ideal for neurodegeneration, autophagy, and protein homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	Ataxin-3 (ATXN3) Rabbit pAb
Remarks/Alias	AT3; JOS; MJD; ATX3; MJD1; SCA3; Ataxin-3 (ATXN3)
Species	Human
Gene ID	4287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-361 of human Ataxin-3 (ATXN3) (NP_004984.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P54252
Protein Weight	41kDa
Shipping	Ice bag