

Anti-Ataxin-3 (Atxn3) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P54252

Item Number: CM0013674



Introduction

Rabbit polyclonal antibody targeting human Ataxin-3 (ATXN3), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying Machado-Joseph disease and ubiquitin pathways.

[Learn More](#)

Parameter	Value
Product Name	Ataxin-3 (ATXN3) Rabbit pAb
Remarks/Alias	AT3; JOS; MJD; ATX3; MJD1; SCA3; Ataxin-3 (ATXN3)
Species	Human
Gene ID	4287
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-361 of human Ataxin-3 (ATXN3) (NP_004984.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P54252
Protein Weight	41kDa
Shipping	Ice bag