


Anti-Atf2 Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, Elisa - P15336

Item Number: CM0022025



Introduction

KO Validated rabbit polyclonal antibody against ATF2 (CREB2, HB16). Suitable for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studying transcriptional regulation, apoptosis, and DNA damage response pathways.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATF2 Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; F2
Species	Human
GeneID (Human)	1386
GeneID	1386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-190 of ATF2 (NP_001871.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag