

Anti-Atf2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - P15336

Item Number: CM0004051



Introduction

High-specificity rabbit monoclonal antibody targeting human ATF2 (CREB-2), validated for WB, IHC-P, IP, ELISA. Cross-reacts with mouse and rat. Suitable for transcriptional regulation and DNA damage response studies.

[Learn More](#)

Parameter	Value
Product Name	ATF2 Rabbit mAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; ATF2
Species	Human
Gene ID (Human)	1386
Gene ID	1386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-340 of human ATF2 (NP_001871.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag