

Anti-Atf2 Rabbit Polyclonal Antibody For Wb, Ihc-P - P15336

Item Number: CM0015577



Introduction

High-quality rabbit polyclonal ATF2 antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat ATF2 protein (~55kDa). Ideal for transcription regulation and DNA damage response research.

[Learn More](#)

Parameter	Value
Product Name	ATF2 Rabbit pAb
Remarks/Alias	HB16; CREB2; TREB7; CREB-2; CRE-BP1; ATF2
Species	Human
GeneID (Human)	1386
GeneID	1386
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human ATF2 (NP_001871.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15336
Protein Weight	55kDa
Shipping	Ice bag