

Anti-Atf3 Rabbit Monoclonal Antibody For Wb - P18847

Item Number: CM0012090

Introduction

ATF3 rabbit monoclonal antibody validated for WB, IP, ChIP, ELISA. Targets human ATF3 (Gene ID 467), cross-reacts with mouse. Recombinant immunogen (aa 1-77). Ideal for transcription and stress research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATF3 Rabbit mAb
Remarks/Alias	ATF3
Species	Human
GeneID (human)	467
GeneID	467
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-77 of human ATF3 (NP_001665.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ChIP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18847
Protein Weight	12-21 kDa
Shipping	Ice bag