

Anti-Atf3 Rabbit Polyclonal Antibody For Wb, Elisa - P18847

Item Number: CM0019411

Introduction

Anti-ATF3 Rabbit Polyclonal Antibody, validated for WB and ELISA, reacts with human, mouse, rat ATF3. Ideal for transcription regulation and stress response studies. P18847.

[Learn More](#)

Parameter	Value
Product Name	ATF3 Rabbit pAb
Remarks/Alias	ATF3
Species	Human
Gene ID (Human)	467
Gene ID	467
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human ATF3 (NP_001665.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18847
Protein Weight	21 kDa
Shipping	Ice bag