

Anti-Atf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P18847

Item Number: CM0017091

Introduction

Rabbit polyclonal antibody specific for ATF3, validated for WB, IHC-P, and ELISA. Suitable for studying stress-induced transcription regulation. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	ATF3 Rabbit pAb
Remarks/Alias	ATF3
Species	Human
Gene ID (human)	467
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human ATF3 (NP_001665.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18847
Protein Weight	21kDa
Shipping	Ice bag