

Anti-Atf4 Polyclonal Antibody For Wb, If/Icc, Elisa - P18848

Item Number: CM0014132



Introduction

High-quality ATF4 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Ideal for integrated stress response, ER stress, and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	ATF4 Rabbit pAb
Remarks/Alias	CREB2; TXREB; CREB-2; TAXREB67; ATF4
Species	Human
GeneID (human)	468
GeneID	468
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P18848
Protein Weight	39kDa
Shipping	Ice bag