

Anti-Atf6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Chip, Elisa - P18850

Item Number: CM0010109

Introduction

High-specificity rabbit monoclonal antibody targeting human ATF6, validated for WB, IHC-P, IP, ChIP, and ELISA. Reacts with human and mouse. Ideal for UPR and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	ATF6 Rabbit mAb
Remarks/Alias	ACHM7; ATF6A; ATF6
Species	Human
Gene ID (Human)	22926
Gene ID	22926
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ChIP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18850
Protein Weight	75 kDa
Shipping	Ice bag