

Anti-Atf6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P18850

Item Number: CM0017834

Introduction

High-quality rabbit polyclonal antibody targeting human ATF6 (ATF6A). Validated for WB, IF/ICC, IHC-P, and ELISA with reactivity across human, mouse, and rat. Ideal for unfolded protein response and ER stress research.

[Learn More](#)

Parameter	Value
Product Name	ATF6 Rabbit pAb
Remarks/Alias	ACHM7; ATF6A; ATF6
Species	Human
GenelD (human)	22926
GenelD	22926
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18850
Protein Weight	75 kDa
Shipping	Ice bag