

Anti-Atf7 Rabbit Polyclonal Antibody For Wb, Elisa - P17544

Item Number: CM0026486

Introduction

Rabbit polyclonal antibody targeting human ATF7, validated for Western blot and ELISA. ATF7 is a stress-responsive chromatin regulator involved in innate immunity, adipocyte differentiation, and telomere maintenance.

[Learn More](#)

Parameter	Value
Product Name	ATF7 Rabbit pAb
Remarks/Alias	ATFA; ATF7
Species	Human
GeneID (Human)	11016
GeneID	11016
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17544
Protein Weight	52kDa
Shipping	Ice bag