


## Anti-Atf7Ip Rabbit Polyclonal Antibody For Wb, Elisa - Q6Vmq6

Item Number: CM0016944



### Introduction

Rabbit polyclonal antibody against human ATF7IP, validated for WB and ELISA. Cross-reacts with mouse. Useful for studying transcriptional regulation, chromatin remodeling, and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATF7IP Rabbit pAb                                                                                                    |
| Remarks/Alias    | AM; MCAF; p621; MCAF1; ATF-IP; ATF7IP1; ATF7IP                                                                       |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 55729                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 525-725 of human ATF7IP (NP_060649.3). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | Q6VMQ6                                                                                                               |
| Protein Weight   | 136kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |