

Anti-Atg10 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9H0Y0

Item Number: CM0008394

Introduction

High-quality rabbit monoclonal antibody against human ATG10, validated for WB, IF/ICC, and ELISA. Reacts with human and mouse ATG10, ideal for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	ATG10 Rabbit mAb
Remarks/Alias	APG10; APG10L; pp12616; ATG10
Species	Human
GeneID (human)	83734
GeneID	83734
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 121-220 of human ATG10 (Q9H0Y0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0Y0
Protein Weight	25kDa
Shipping	Ice bag