

Anti-Atg10 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H0Y0

Item Number: CM0027322

Introduction

Rabbit polyclonal antibody targeting ATG10, an E2-like enzyme that catalyzes ATG12-ATG5 conjugation during autophagy. Suitable for WB, IHC-P, and ELISA in human, mouse, and rat samples. Immunogen: synthetic peptide within amino acids 120-220.

[Learn More](#)

Parameter	Value
Product Name	ATG10 Rabbit pAb
Remarks/Alias	APG10; APG10L; pp12616; ATG10
Species	Human
GeneID (Human)	83734
GeneID	83734
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 120-220 of human ATG10 (NP_113670.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0Y0
Protein Weight	25kDa
Shipping	Ice bag