

Anti-Atg12 Polyclonal Antibody For Wb, If/Icc, Elisa - 094817

Item Number: CM0014212



Introduction

High-quality rabbit polyclonal antibody against human ATG12 (autophagy-related protein 12). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for autophagy and host-virus interaction research.

[Learn More](#)

Parameter	Value
Product Name	ATG12 Rabbit pAb
Remarks/Alias	APG12; FBR93; APG12L; HAPG12; ATG12
Species	Human
GeneID	9140
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATG12 (NP_004698.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94817
Protein Weight	15kDa
Shipping	Ice bag