

## Anti-Atg12 Rabbit Monoclonal Antibody For Wb, Elisa - 094817

Item Number: CM0002160



### Introduction

Rabbit monoclonal antibody against human ATG12 (APG12L) for Western blot and ELISA. Validated for autophagy research with cross-reactivity to human and mouse samples.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | ATG12 Rabbit mAb                                                                                                     |
| Remarks/Alias    | APG12; FBR93; APG12L; HAPG12; ATG12                                                                                  |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 9140                                                                                                                 |
| Gene ID          | 9140                                                                                                                 |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 41-140 of human Apg12 (094817). |
| Source           | Rabbit                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | 094817                                                                                                               |
| Protein Weight   | 15 kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                              |