

Anti-Atg12 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - O94817

Item Number: CM0007028



Introduction

This ATG12 rabbit monoclonal antibody detects human and mouse ATG12 protein (15 kDa) in WB, IHC-P, IP, ELISA. Ideal for autophagy and host-virus interaction studies.

[Learn More](#)

Parameter	Value
Product Name	ATG12 Rabbit mAb
Remarks/Alias	APG12; FBR93; APG12L; HAPG12; ATG12
Species	Human
GeneID (Human)	9140
GeneID	9140
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ATG12 (NP_004698.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	O94817
Protein Weight	15kDa
Shipping	Ice bag