


Anti-Atg13 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O75143

Item Number: CM0017636



Introduction

Rabbit polyclonal antibody against ATG13 (UniProt O75143), suitable for Western blot, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat ATG13. Supports autophagy and mitophagy studies.

[Learn More](#)

Parameter	Details
Product Name	ATG13 Rabbit pAb
Remarks/Alias	KIAA0652; PARATARG8; ATG13
Species	Human
Gene ID (Human)	9776
Gene ID	9776
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-480 of human ATG13 (NP_001136145.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75143
Protein Weight	57 kDa
Shipping	Ice bag