

Anti-Atg14 Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Zne5

Item Number: CM0028676



Introduction

Human ATG14 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, rat. Detects 55 kDa Beclin1-associated autophagy regulator, ideal for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	ATG14 Rabbit pAb
Remarks/Alias	ATG14L; BARKOR; KIAA0831; ATG14
Species	Human
GeneID (Human)	22863
GeneID	22863
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-440 of human ATG14 (NP_055739.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6ZNE5
Protein Weight	55kDa
Shipping	Ice bag