

Anti-Atg16L1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q676U5

Item Number: CM0013560



Introduction

Rabbit polyclonal antibody targeting ATG16L1, key for autophagy. Suitable for WB, IHC-P, IF/ICC, ELISA across human, mouse, rat. Ideal for autophagy and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	ATG16L1 Rabbit pAb
Remarks/Alias	IBD10; WDR30; APG16L; ATG16A; ATG16L; ATG16L1
Species	Human
GeneID	55054
Immunogen	A synthetic peptide within amino acids 200-300 of human ATG16L1 (NP_110430.5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q676U5
Protein Weight	68kDa
Shipping	Ice bag