

Anti-Atg16L1 Polyclonal Antibody - Q676U5

Item Number: CM0021794



Introduction

Rabbit polyclonal antibody targeting human ATG16L1 (SWISS: Q676U5). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Ideal for autophagy and inflammatory pathway studies.

[Learn More](#)

Parameter	Value
Product Name	ATG16L1 Rabbit pAb
Remarks/Alias	IBD10; WDR30; APG16L; ATG16A; ATG16L; ATG16L1
Species	Human
GeneID	55054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-300 of human ATG16L1 (NP_110430.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q676U5
Protein Weight	68kDa
Shipping	Ice bag