

Anti-Atg16L1 Rabbit Monoclonal Antibody For Wb, Elisa - Q676U5

Item Number: CM0009550



Introduction

Recombinant rabbit monoclonal antibody targeting human ATG16L1, validated for Western blot and ELISA applications. Reacts with human and mouse species. Ideal for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] ATG16L1 Rabbit mAb
Remarks/Alias	IBD10; WDR30; APG16L; ATG16A; ATG16L; ATG16L1
Species	Human
GeneID (Human)	55054
GeneID	55054
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 6-105 of human ATG16L1 (Q676U5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q676U5
Protein Weight	68 kDa
Shipping	Ice bag