

## Anti-Atg2A Polyclonal Antibody For Wb, Elisa - Q2Taz0

Item Number: CM0026572



### Introduction

Rabbit polyclonal antibody targeting human ATG2A (BLTP4A), suited for Western blot and ELISA. A key lipid transfer protein in autophagosome biogenesis and lipid droplet regulation.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | ATG2A Rabbit pAb                                                                                      |
| Remarks/Alias    | BLTP4A; ATG2A                                                                                         |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 23130                                                                                                 |
| GeneID           | 23130                                                                                                 |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 700-800 of human ATG2A (NP_055919.2) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human                                                                                                 |
| SWISS            | Q2TAZ0                                                                                                |
| Protein Weight   | 213kDa                                                                                                |
| Shipping         | Ice bag                                                                                               |