


Anti-Atg3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nt62

Item Number: CM0002041



Introduction

Anti-ATG3 rabbit monoclonal antibody validated for Western blot, IHC-P, and ELISA detection of human, mouse, and rat ATG3. Supports reliable autophagy research with high specificity.

[Learn More](#)

Parameter	Specification
Product Name	ATG3 Rabbit mAb
Remarks/Alias	APG3; APG3L; hApg3; PC3-96; APG3-LIKE
Species	Human
GeneID (Human)	64422
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Apg3 (Atg3) (Q9NT62)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NT62
Protein Weight	40kDa
Shipping	Ice bag