


Anti-Atg3 Rabbit Polyclonal Antibody - Q9Nt62

Item Number: CM0029201



Introduction

High-quality rabbit polyclonal antibody against human ATG3, validated for WB, IF/ICC, and ELISA. Suitable for autophagy research with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	ATG3 Rabbit pAb
Remarks/Alias	APG3; APG3L; hApg3; PC3-96; APG3-LIKE; ATG3
Species	Human
GeneID (human)	64422
GeneID	64422
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NT62
Protein Weight	36kDa
Shipping	Ice bag