


Anti-Atg4A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wyn0

Item Number: CM0020530



Introduction

Rabbit polyclonal antibody targeting human cysteine protease ATG4A (APG4A/AUTL2). Suitable for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Raised against recombinant protein (aa 210-398). For autophagy and delipidation studies.

[Learn More](#)

Parameter	Value
Product Name	ATG4A Rabbit pAb
Remarks/Alias	APG4A; AUTL2; HsAPG4A; ATG4A
Species	Human
GeneID (human)	115201
GeneID	115201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-398 of human ATG4A (NP_443168.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WYN0
Protein Weight	45kDa
Shipping	Ice bag