


Anti-Atg4B Polyclonal Antibody For Wb, Elisa - Q9Y4P1

Item Number: CM0028791



Introduction

Rabbit polyclonal antibody targeting human ATG4B (APG4B, AUTL1), a key autophagy protease. Validated for Western Blot and ELISA. Suitable for studying autophagy, mitophagy, and ATG8 conjugation pathways.

[Learn More](#)

Parameter	Value
Product Name	ATG4B Rabbit pAb
Remarks/Alias	APG4B; AUTL1; HsAPG4B; ATG4B
Species	Human
Gene ID (Human)	23192
Gene ID	23192
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human ATG4B (NP_037457.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y4P1
Protein Weight	44kDa
Shipping	Ice bag