


Anti-Atg4B Rabbit Monoclonal Antibody For Wb, If, Elisa - Q9Y4P1

Item Number: CM0004589



Introduction

Rabbit monoclonal antibody targeting human ATG4B (Q9Y4P1), validated for WB, IF/ICC, IF-P, and ELISA. Cross-reacts with mouse and rat. Essential tool for studying autophagy, mitophagy, and ATG8 processing.

[Learn More](#)

Parameter	Value
Product Name	ATG4B Rabbit mAb
Remarks/Alias	APG4B; AUTL1; HsAPG4B; ATG4B
Species	Human
GeneID (Human)	23192
GeneID	23192
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ATG4B (Q9Y4P1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4P1
Protein Weight	44kDa
Shipping	Ice bag