

Anti-Atg4C Polyclonal Antibody For Wb, Elisa - Q96Dt6

Item Number: CM0027946

Introduction

Rabbit polyclonal antibody targeting human ATG4C, validated for WB and ELISA. ATG4C is a cysteine protease essential for autophagy, processing ATG8 proteins. Ideal for autophagy research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	ATG4C Rabbit pAb
Remarks/Alias	APG4C; AUTL1; AUTL3; APG4-C; HsAPG4C; ATG4C
Species	Human
GeneID (Human)	84938
GeneID	84938
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ATG4C (NP_116241.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96DT6
Protein Weight	52kDa
Shipping	Ice bag