

Anti-Atg4D Polyclonal Antibody For Wb And Elisa - Q86TL0

Item Number: CM0030284

Introduction

Anti-ATG4D polyclonal antibody validated for WB and ELISA, reactive with human and rat ATG4D. Cysteine protease involved in autophagy; suitable for investigating autophagy and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	ATG4D Rabbit pAb
Remarks/Alias	APG4D; AUTL4; APG4-D; HsAPG4D; ATG4D
Species	Human
Gene ID (Human)	84971
Gene ID	84971
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 265-474 of human ATG4D (NP_116274.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q86TL0
Protein Weight	53kDa
Shipping	Ice bag