

Anti-Atg4D Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q86TL0

Item Number: CM0021021

Introduction

Rabbit polyclonal antibody targeting human ATG4D (APG4D/AUTL4), validated for WB, IF/ICC, and ELISA. Reacts with human and mouse. ATG4D is a cysteine protease essential for autophagy and mitophagy.

[Learn More](#)

Parameter	Value
Product Name	ATG4D Rabbit pAb
Remarks/Alias	APG4D; AUTL4; APG4-D; HsAPG4D; ATG4D
Species	Human
GeneID (Human)	84971
GeneID	84971
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human ATG4D (NP_116274.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86TL0
Protein Weight	53kDa
Shipping	Ice bag