

Anti-Atg5 Monoclonal Antibody For Wb, If-P, Ihc-P, Ip, Elisa - Q9H1Y0

Item Number: CM0002087



Introduction

KO validated anti-ATG5 rabbit monoclonal antibody for WB, IF-P, IHC-P, IP, ELISA. Species reactivity: human, mouse, rat. Ideal for autophagy, apoptosis, and lymphocyte development studies.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] ATG5 Rabbit mAb
Remarks/Alias	ASP; APG5; APG5L; hAPG5; SCAR25; APG5-LIKE; G5
Species	Human
Gene ID (Human)	9474
Gene ID	9474
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1Y0
Protein Weight	32kDa
Shipping	Ice bag