

Anti-Atg5 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9H1Y0

Item Number: CM0017101



Introduction

High-quality rabbit polyclonal antibody recognizing human ATG5 (Autophagy protein 5) with cross-reactivity to mouse and rat. Ideal for Western blot, IHC-P, and ELISA applications in autophagy research.

[Learn More](#)

Parameter	Value
Product Name	ATG5 Rabbit pAb
Remarks/Alias	ASP; APG5; APG5L; hAPG5; SCAR25; APG5-LIKE; ATG5
Species	Human
Gene ID (Human)	9474
Gene ID	9474
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-275 of human ATG5 (NP_004840.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1Y0
Protein Weight	32kDa
Shipping	Ice bag