

Anti-Atg7 Polyclonal Antibody For Ihc-P And Elisa - O95352

Item Number: CM0025806

Introduction

Rabbit polyclonal antibody targeting human ATG7 (O95352), validated for IHC-P and ELISA. Ideal for autophagy research, mitophagy, and biomarker detection in human samples.

[Learn More](#)

Parameter	Details
Product Name	ATG7 Rabbit pAb
Remarks/Alias	GSA7, APG7L, SCAR31, APG7-LIKE, ATG7
Species	Human
GeneID (Human)	10533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-676 of human ATG7 (NP_001129503.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	O95352
Protein Weight	78kDa
Shipping	Ice bag