

Anti-Atg7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95352

Item Number: CM0015385

Introduction

High-quality rabbit polyclonal antibody against human ATG7 (autophagy-related 7). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	ATG7 Rabbit pAb
Remarks/Alias	GSA7; APG7L; SCAR31; APG7-LIKE; ATG7
Species	Human
Gene ID (Human)	10533
Gene ID	10533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-676 of human ATG7 (NP_006386.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95352
Protein Weight	78 kDa
Shipping	Ice bag