

Anti-Atg8 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Xy24

Item Number: CM0021783

Introduction

Rabbit polyclonal antibody against wheat ATG8, validated for Western blot and ELISA. Reacts with Arabidopsis, maize, and wheat. Suitable for autophagy research in plants.

[Learn More](#)

Parameter	Value
Product Name	ATG8 Rabbit pAb
Species	Triticum aestivum
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of Triticum aestivum ATG8 (Q7XY24)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Arabidopsis thaliana, Zea mays, Triticum aestivum
SWISS	Q7XY24
Protein Weight	14kDa
Shipping	Ice bag