


Anti-Atg9A Polyclonal Antibody For Wb - Q7Z3C6

Item Number: CM0026812



Introduction

High-quality rabbit polyclonal antibody against ATG9A, raised against a recombinant immunogen (aa 180-290). Validated for WB and ELISA. Detects human, mouse, and rat ATG9A. Ideal for autophagy research.

[Learn More](#)

Parameter	Specification
Product Name	ATG9A Rabbit pAb
Remarks/Alias	mATG9, APG9L1, MGD3208, ATG9A
Species	Human
Gene ID (Human)	79065
Gene ID	79065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human ATG9A (NP_076990.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z3C6
Protein Weight	94kDa
Shipping	Ice bag