


Anti-Atg9B Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q674R7

Item Number: CM0027013



Introduction

Rabbit polyclonal antibody against human ATG9B, validated for Western blot, immunofluorescence, and ELISA. Suitable for human, mouse, and rat samples. Ideal for autophagy and cell death research.

[Learn More](#)

Parameter	Value
Product Name	ATG9B Rabbit pAb
Remarks/Alias	SONE; APG9L2; NOS3AS; ATG9B
Species	Human
Gene ID	285973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 715-924 of human ATG9B (NP_001303985.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q674R7
Protein Weight	101kDa
Shipping	Ice bag