

Anti-Atic Rabbit Polyclonal Antibody For Wb, Elisa - P31939

Item Number: CM0025225



Introduction

Rabbit polyclonal antibody targeting human ATIC, a bifunctional purine biosynthesis enzyme. Validated for Western blot and ELISA applications. Ideal for purine metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ATIC Rabbit pAb
Remarks/Alias	PURH; AICAR; AICARFT; IMPCHASE; HEL-S-70p; ATIC
Species	Human
GeneID (Human)	471
GeneID	471
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human ATIC (NP_004035.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P31939
Protein Weight	65kDa
Shipping	Ice bag