

Anti-Atl2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nhh9

Item Number: CM0019167

Introduction

ATL2 rabbit polyclonal antibody, cross-reacts with human, mouse, and rat, suitable for WB, IF/ICC, and ELISA. Useful for investigating endoplasmic reticulum network maintenance and GTP-dependent membrane fusion.

[Learn More](#)

Parameter	Value
Product Name	ATL2 Rabbit pAb
Remarks/Alias	aip-2; ARL3IP2; ARL6IP2; atlastin2; ATL2
Species	Human
GeneID	64225
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 390-470 of human ATL2 (NP_001129145.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NHH9
Protein Weight	66kDa
Shipping	Ice bag