


Anti-AtI3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Dd88

Item Number: CM0015628



Introduction

Rabbit polyclonal antibody against human ATL3 (Atlastin-3), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for ER dynamics and membrane fusion research.

[Learn More](#)

Parameter	Value
Product Name	ATL3 Rabbit pAb
Remarks/Alias	HSN1F; ATL3
Species	Human
Gene ID (Human)	25923
Gene ID	25923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 360-440 of human ATL3 (NP_056274.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6DD88
Protein Weight	61kDa
Shipping	Ice bag